

Integration Command for AD131

---oOo---

1. Introduction

This command allows a user to set the integration time in microseconds (μs) instead of using two gain commands.

- L-Set the integration variable (gain)
< L > (1 - 255) < 1 - 255 >
- X-Extension to gain range. [1-255] (x-gain)
< X > (present extended gain value) < new extended gain value: 1 - 255 >

This command is used the same purpose as the gain and x-gain commands that are to increase or decrease the exposure time of the detector.

2. Command Structure

To AD131: < B >

From AD131: (High Byte)(Middle Byte)(Low Byte)

To AD131: <High Byte><Middle Byte><Low Byte>

This command is in effect when the integration time sets to $88\mu\text{s}$ or longer time.

Otherwise, the gain and x-gain will be in effect. The default integration is $88\mu\text{s}$.

- Example 1:
 - i. if a user sets the integration to $2005\mu\text{s}$, then the integration time is $2005\mu\text{s}$.
 - ii. If a user sets the integration to $23\mu\text{s}$, the default gain and x-gain is in effect.

- Example 2:

The current integration is $23\mu\text{s}$, and a user wants to set the integration to $2005\mu\text{s}$.

Then the user does the followings:

To AD131: B

From AD131: 0, 0, 23

To AD131: 0, 7, 213

Notes:

- $2005_{10} = 7D5_{16} = 00\ 07\ D5_{16}$
- $07_{16} = 7_{10}$
- $D5_{16} = 213_{10}$